

## IPTG Induction of recombinant protein expression in Bacteria

### Protocol:

- 1) At the end of the day before the lab day, pick a colony and grow a 5ml overnight culture at 30/37 °C in LB plus the right antibiotic in a snap cap test tube, in a rotator or shaker. Note: This step can be performed by the TAs if none of your group member is available to do it.
- 2) In the next morning, check the absorbance of the culture at 600 nm. Inoculate a new 5 ml culture with the overnight culture. Use a dilution factor that will make the new culture start from  $OD_{600}=0.15$ . It will help if you incubate the medium in a test tube for the new culture in the same incubator so that it will be at the right temperature when you inoculate it.
- 3) Start checking the OD of the new culture 2 hours after the inoculation. When the  $OD_{600}$  reaches 0.6 (it would be fine up the  $OD_{600}=1.2$ ), remove 0.5 ml culture from the tube and place it 1.7ml tubes. Spin the tube at 5000 rpm for 1 min. Remove the supernatant and add 100  $\mu$ l a 2X SDS sample buffer. Boil it for 5 min. and freeze it for later use. THIS IS THE UNINDUCED CONTROL.
- 4) Add 4  $\mu$ l of a 1000x IPTG stock to the culture to start the induction of the protein expression.
- 5) After 3 hours (it would be fine up to 7 hours), transfer 0,5 ml of the culture to a 1.7 ml tubes and follow the same procedure in step 3 to prepare the INDUCED SAMPLE. Store the samples in the -20C° freezer for the next lab session. Note: for our class this may mean that some of your group members will need to come and harvest the cells at the end of the lab day.

### Lab day 2. Electrophoresis:

- 6) Follow the protocols in another handout to prepare an SDS-gel.
- 7) Load your samples (the TAs will show you how) in the wells and run it